

## Product details

### The feature of tie-on face mask

- 1) Fluid-resistant, high breath ability and splash protection.
- 2) Type IIR quality standard, conforms to EN 14683, CE standard
- 3) The bacterial filtration efficiency (BFE) of 99 percent
- 4) It's pleated design and vertical bands afford excellent facial fit
- 5) Made of a special polypropylene material, non cytotoxic
- 6) Pliable nose shaping bar, and full chin protection.

### Packaging and dimensions of hygiene face mask

50 pcs/box, box size: 18\*10\*10cm

2000 pcs/ctn, ctn size: 52\*38\*42cm

Or as your requirements.

### The trade terms of mouth cover mask

MOQ	100,000 pcs
Price	Negotiable
Payment terms	T/T, L/C, D/A, D/P, Western Union, ESCROW, Paypal
Supply ability	1,000,000 pcs/day, 1*40 HQ container can be finished in 15 days
Delivery time	Normally 30 days after confirmed the order
Sample	Sample is free, can be prepared well in 2 days

## Product display



### The regular shipment way



## FAQ

Q1: What is your output capacity for your major range of the product?

A: Face mask capacity: 1000,000 pieces / day.

Coverall, Surgical gowns, lab coat, bed sheet capacity: 10,000 pieces/day.

Caps capacity: 200,000 pieces/day. Each 1\*40HQ can be finished within 1 week.

Q2: what is your MOQ of each products?

A: Face mask/caps/shoe covers/arm sleeves each MOQ: 50,000 pieces.

Clothes and sheets MOQ:5000 pieces