

1. PRODUCT SIZE:

S M L XL XXL

2. MATERIAL:

PP 30-70GSM

SMS 30-70GSM

SF 50-70CM

PP+PE 30-60GSM

3. MATERIAL TREATMENT:

Material TREATMENT RESULT

ANTI-WATER: YES

ANTI-DUST: YES

ANTI-STATIC: YES/NO

ANTI-BLOOD: YES/NO

ANTI-ALCOHOL: YES/NO

4. COLOR: WHITE, BLUE, DARK BLUE, ORANGE AND SO ON**5. USAGE:**

A. KIMONO be used in the spa and sauna to convenience the body against the water and the dust.

B. KIMONO can be used in the massage to relax your body.

C. KIMONO can be used in the cleaning room, for example the GMP.

D. KIMONO can be used in the hospital.

6. PACKING:

10pc/bag , 100pcs/ctn

7. OEM SERVICE:

A. Can print the band for our customer on the PE bag

B. Can print the band for our customer on the Carton box

8. SAMPLES:

Samples can be sent out within 3 days after you confirmed.

9. DELIVERY DATE:

It can be sent within 15 days after received the deposit.

