

Description of [Disposable Isolation Gowns](#):

1. Directly manufacturer, any colors.
2. ISO13485, CE, PP, PP/PE, SMS, SMMS
3. Capacity: 10,000 pcs/day, 1*40HQ can be finished within 10-15 days;
4. Good barrier to fluid and breathable, OEM packing design, service in 24 hours.

Feature of [Disposable Isolation Gowns](#):

1. Disposable Isolation Gowns is made from high-performance isolation protective material.
2. Widely used in medical safety protective, protective isolation, farm animal husbandry biohazard and soon.
3. With two kinds processes: sewing and ultrasound thermal bonding.

Specification of [Disposable Isolation Gowns](#):

Name	Wholesaler Personal Touch Universal Size Blue Disposable Isolation Gowns in China
Material	PP, PP/PE, SMS, SMMS, CPE
Style	Tie/ loop fastener on collar, tie on waist, elastic/knitted cuff
Size	S-6XL
Color	white, green, blue, yellow, pink, red, grey, black etc.
Packing	Non-sterilized : 5 Pcs/box, 10 bags/ctn EO-sterilized : 1 Pcs/paper plastic bag, 50 PCs/Ctn. Or as your requirement.
Certificate	CE, ISO13485, FDA
MOQ	500 Pcs
Sample	Can be supplied free for your quality checking within 3 days
Payment way	T/T, PayPal, Western union, L/C, D/A, D/P, ESCROW,
Delivery time	Capacity: 10,000 pcs/day, 1*40HQ can be finished within 10-15 days;
FOB port	Wuhan/Shanghai/Ningbo/Guangzhou and as you required

Detail Size of Disposable Isolation Gowns:

Size	S	M	L	XL
Width (cm)	110	115	120	125
Height (cm)	130	137	140	150

Details Pictures of Disposable Isolation Gowns:



