

Product Details:

1. Material: Non woven polypropylene 100%
2. Size:80*190+20cm,90*200+20cm,160*190+20cm,160*200+20cm,180*200+20cm or customized size
3. Elastic: 5mm along side width+30cm on length or customized size
4. Colour: Blue; White ; other size available
5. Packing: 10pcs/bag, 100pcs/carton

Indications

- 1, Widely used for SPA, Hotel, Massage, Salon, Medical one time use
- 2, To prevent the spread of body fluids and blood splash content, and isolate the dust, particle, alcohol, blood, bacteria and virus invading.

Specifications:

Material	Non woven polypropylene 100%
Color	White, blue
Style	elastics on width
Size	80*190+20cm,90*200+20cm,160*190+20cm,160*200+20cm,180*200+20cm or customized size
packing	10 pcs/box, 100pcs/ctn
Weight	25-65gsm or customized
Package	PE bag coated, then paper carton package or customized
Feature	Durable, soft, disposable
Usage	SPA, Hotel, Massage, Salon, Medicinal, Hospital, Home

Our Services:

1. Provide the highest quality products

Professional QC team strictly control each step of the production until loading.

2. The fastest delivery time

If urgent, 1*40HQ can be finished within 5-7days!

Normally, 1*40HQ can be finished within 10-15 days!

Busy season, 1*40HQ can be finished within 25-30 days!

3. The most reasonable price

Continuously invest in improving our skilled workers, products and processes to minimize the cost of each step!

4.The sweet and quick service

Our workers will be always so patiently making a response of your any questions within 5-10 minutes.

5.We can also help to source connected products.

Besides all kinds of disposable medical products, we can also help your source the connected products as your demand!

<http://https://www.china-disposable.com/products/CE-ISO-approved-liquid-proof-one-time-use-non-woven-bed-cover-with-elastic-ends.html#.YHAKNOgzaUk>