

Product details

best quality supplier for disposable anti-static lab coat

SMS medical disposable clothes made from hydrophobic SMS material ,Latex-free; abrasion-resistant; low lint; with a high level of fluid repellency; a good barrier to blood, body fluids and pathogens.

Available color: light blue, blue, dark blue

Specification of best quality supplier for disposable anti-static lab coat:

Product:	best quality supplier for disposable anti-static lab coat,SMS medical disposable clothes , wholesales Disposable hospital SMS Clothing
Material:	SMS
Grammage:	50G
Style:	long sleeves; knit cuff ; knitted collar with three pockets;with snap button closure
Application:	for medical, pharmaceutical, laboratorial, workshop sector
Color:	Blue,
MOQ	5000pcs
Packaging:	5 or 10pcs/bag; 50 or 100pcs/ctn;
Storage:	Stored in dry, humidity below 80%, ventilated, non-corrosive gases warehouse
Terms of payment:	T/T, PayPal, Western union, L/C,D/A,D/P,ESCROW,
Other service:	OEM and Customized Logo/brand
Certification:	CE, ISO13485
Delivery time	Capacity: 10,000pcs/day, 1*40HQ can be finished within 10-15 days;

Product display



orient Honest
MEDICAL INTERNATIONAL






orient Honest
MEDICAL INTERNATIONAL



Our business process



The regular shipment way



FAQ

Q: 1.Are you a manufacturer or a trading company, and when was your company established?

A: Yes, we are directly manufacturer, our factory was established in 2003, we have our own exported rights company in 2013.

Q: 2.What is the annual turnover for your company?

A: Normally it is around USD18 million-20 million, it is still in the growing.

Q: 3. Do your company have any quality management system such as ISO, EN or CE certificate(s)?

A: Yes, we have CE,ISO13485,FDA. REACH and NELSON report.

Q: 4.Can I get the free samples for quality checking? How long time can I get it?

A: Yes, our samples can be supplied free within 3-7 days.