

Colored lab coats with knitted cuff:

colored lab coats with knitted cuff, factory custom made pink lab coats, good quality vistor coat made from PP+PE material, Latex-free; abrasion-resistant; completely impervious fluid and alcohol repellent.

Available color: pink, white, green, blue, yellow, etc

Material weight: 25-65gsm

Colored lab coats with knitted cuff size chart:

Size	S	M	L	XL	2XL	3XL
Height (CM)	103	111	115	120	124	128
chest circumference(CM)	116	126	132	34	144	154
Sleeve length(CM)	57	57	57	58	58	58


Specification of colored lab coats with knitted cuff:

Name	colored lab coats with knitted cuff, factory custom made pink lab coats, good quality vistor coat manufacturers
Material	PP+PE
Grammage	25-65gsm/ 30-65gsm/45-65gsm or as your request
Style	With/without pocket, single/double collar, Velcro/zipper/button on the front ,elastic/knitted cuff
Size	104x148cm,108x152cm,112x156cm,115x160cm or as your request
Color	white, green, blue,
packing	5pcs/bag, 10bags/carton or as your request
packing design	All inner bag and carton printing can be as your requirement
Sample	Can be supplied free for your quality checking within 3 days
Certifications	CE / ISO13485 / FDA

Our Advantage:

1. Samples and customization available.
2. Professional QC to control the product quality. All the products comply with international quality standard such as CE,ISO13485.
3. Have enough experience and gained a high reputation from our customers all the world. Such as America, South America, Europe, Southeast Asia, The Middle East etc.
4. Quick payment, prompt delivery and competitive price.
5. Stable quality and best after-sale service.
6. OEM service is available.

Pictures of colored lab coats with knitted cuff:





orient Honest
MEDICAL INTERNATIONAL