

## Product details

### Introduction of custom anti-static disposable gown

Custom anti-static disposable gown ,China the best disposable doctor gown,China factory disposable surgeon gowns is made from hydrophobic polypropylene material with PE coated, or Spunlace material with PE coated.With the good function for anti alcohol,anti blood,waterproof,etc

### The packaging of custom anti-static disposable gown

Sterile: 1 pcs/paper plastic bag,50 pcs/Ctn

Non-sterilize: 5 pcs/bag, 50 pcs/ctn,Or as your requirements.

### Specifications of custom anti-static disposable gown

Name:	custom anti-static disposable gown ,China the best disposable doctor gown,China factory disposable surgeon gowns
Material:	PP+PE45G
Grammage:	45 gsm
Closure:	with belt closure, or with Velcro closing
Application:	for medical, surgical, pharmaceutical, laboratorial, industrial sector
Size:	S/M/L/XL/XXL/XXXL,
Color:	Green, etc, optional
Packing:	Normal packing:1pcs/polybag,50pcs/carton
Storage:	Stored in dry, humidity below 80%, ventilated, non-corrosive gases warehouse
Certification:	CE, ISO13485,FDA

## Product display




orient Honest  
MEDICAL INTERNATIONAL







## Our business process



## The regular shipment way



## FAQ

Q: 1. Are you a manufacturer or a trading company, and when was your company established?

A: Yes, we are directly manufacturer, our factory was established in 2003, we have our own exported rights company in 2013.

Q: 2. What is the annual turnover for your company?

A: Normally it is around USD18 million-20 million, it is still in the growing.

Q: 3. Do your company have any quality management system such as ISO, EN or CE certificate(s)?

A: Yes, we have CE, ISO13485, FDA, REACH and NELSON report.

Q: 4. Can I get the free samples for quality checking? How long time can I get it?

A: Yes, our samples can be supplied free within 3-7 days.