

Disposable PP SMS PP+PE gown level1 level2 level3 protective blue gown

Disposable gown is made from high-performance isolation protective material.

Widely used in medical safety protective, aseptic workshop, protective isolation, mining electronics factory, food factory farm animal husbandry biohazard and so on.

With two kinds processes: sewing and ultrasound thermal bonding.

Material features :

1. PP disposable gown

It is made from hydrophobic polypropylene material, Latex-free; abrasion-resistant; low lint; with a high level of fluid repellency.

Color: white, green, blue, yellow, orange and so on

Material weight: 20-65gsm.

2. PP+PE disposable gown

It is made from PP+PE material, Latex-free; abrasion-resistant; completely impervious fluid and alcohol repellent.

Color: white, green, blue, yellow.

Material weight: 40-65gsm.

3. SMS disposable gown (Standard and Reinforced):

It is made from hydrophobic SMS/PP/PP+PE material, Latex-free; abrasion-resistant; low lint; with a high level of fluid repellency; a good barrier to blood, body fluids and pathogens.

With standard and reinforced style. Reinforced style with extra protection reinforcement on sleeve & chest which can be completely impervious fluid and alcohol repellent.

Available color: white, blue.

Material weight: 35-65gsm.

Disposable gown style:

Tie/ hook and loop fastener on collar, tie on waist, elastic/knitted cuff

Disposable gown Specifications

Description	Disposable sterile reinforced surgical gowns with reasonable price
Material	PP, PP+PE, SMS, CPE

Weight	PP 20-65gsm; PP+PE 40-65 gsm; SMS 35-65 gsm (comfortable,without glass fibres,latex free.)
Style	tie/ hook and loop fastener on collar,tie on waist,elastic/knitted cuff
Size	120x130cm,130x150cm,150x160cm,150x190cm, 170x170cm
Color	white,green,blue,yellow
packing	non sterile,10 pcs/bag,10 bags/ctn. sterile,1 pcs/sterile bag,100 pcs/ctn
packing design	all inner box and carton printing can be as your design
carton size	58X28X48 cm or as your requirement
Applications	hospital clinical medical personnel and patients. dust-free workshop,laboratory, food industry, electronic manufactures and so on
MOQ	10000 pcs
Sample	can be supplied free for your quality checking within 3 days
Payment term	T/T,PayPal
Delivery time	as your quantity(we can finish 20,000 pcs per week)
FOB port	Wuhan/Shanghai/Ningbo/Guangzhou