

Disposable PP work clothes non woven one time use pink colour lab [coat](#) is made from hydrophobic SMS, PP, PP+PE. CPE, Microporous material, Latex-free; Abrasion-resistant; Low lint;

Available color: Pink, white, green, blue, yellow, red, black etc. Material weight: 40-65GSM

Lab coats size: S 96*140*57cm M 100*144*57 cm L 104*148*57 XL 108*152*58

Disposable lab coats products Advantages

1. With a high level of fluid repellency.
2. Latex-free; Abrasion-resistant; Completely impervious fluid. A good barrier to blood, body fluids and pathogens.
3. With standard and reinforced style. Reinforced style with extra protection reinforcement on sleeve & chest which can be completely impervious fluid and alcohol repellent.
4. Designed for optimal protection, comfortable wearing.
5. Capacity: 5000pcs/day, 1*40HQ can be finished within 10-15 days;

Disposable lab coats indications

Disposable lab coats are widely used in medical safety protective gowns, aseptic workshop gowns, protective isolation,

Mining electronics factory, food factory farm animal husbandry biohazard and so on.

Disposable lab coats normal packing: 5 pcs/bag, 10 bags/ctn;

Name	disposable PP work clothes non woven one time use pink colour lab coat
Material	PP, PP/PE, SMS, SMMS, Microporous
Style	Knitted collar, knitted cuff , 5 buttons , 3pockets
Size	S-5XL
Color	Pink, white, green, blue, yellow, red, black etc.
packing	5pcs/box, 10bags/ctn □ 10pcs/box, 10bags/carton, or as customer's requirement
packing design	all inner box and carton printing can be as your design
MOQ	5000 pcs
Sample	Lab coats sample can be supplied free for your quality checking within 3 days
Payment way	T/T, PayPal, Western union, L/C, D/A, D/P, ESCROW,

 orient Honest
MEDICAL INTERNATIONAL



