

Description of disposable child lab coat:

1. Name: disposable child lab coat
2. Application: Hospital, hotel, dustproof place, food/electronic/chemical workshop
3. Material:pp 45g
4. Color: White,etc
5. Style: knitted cuff/collar, front button
6. Size: S-XXXL sizes available
7. Packing: 5PCS/polybag, 50PCS/carton and as you require
- 8.1*20" FCL: 50000PCS (tolerance as different sizes and materials)
9. Mini quantity: 10, 000PCS
10. Certificate: CE, FDA, ISO13485

Features of disposable child lab coat:

- 1.Non-irritating,antibacterial property ,good tensile strength,Lint-free,Superior breathability and anti-static,impact resistance,soft,comfortable.
- 2.Against splatter of chemical material and solid dust.
- 3.Usually nonwoven lab coat with the effect of dust-free and anti-virus be widely used public place,medical industry,environment cleaning,electronic factory,food processing,school and hospital.

Our advantage:

1. Direct manufacturer, factory price ;
2. 10,000 Pcs/day,we can finish 1*40 HQ container in about 15days;
3. CE ISO FDA approved , stable quality
4. Service in 24 hours, sample can be finished in 2 days, sample free.

Specifications of disposable child lab coat

Material	PP nonwoven
Fabric weight	45g,or Could as your request
Style	Knitted collar, knitted cuff, without three pockets
Size	S-6XL
Color	white,green,blue,yellow,pink,red, gray,black etc.
Packing	5pcs/box,50 bags/ctn, or as your requirement
Packing design	all inner bag and carton printing can be as your design

Applications	dust free workshop, hospital, laboratory, food industry, electronic manufactures and so on
--------------	--

Pictures of disposable child lab coat:





orient Honest
MEDICAL INTERNATIONAL

