

## Description of disposable clip cap:

Material : PP nonwoven

Color: white,green,blue,pink,yellow,pink,orange,purple,red (all colors can be available)

Style: single elastic or double elastic

Weight: 10g-30g

Size: 18",19",20",21",28"

Packing details: 100 pcs/bag,20 pcs/ctn

Certification: CE / ISO13485 / FDA

Place of Origin: China

## Specification of disposable clip cap:

Material	PP non woven
Weight	10-30g
Style	Elastic (latex free), single and double elastic style
Size	18",19",21",24",28"
Color	white,green,blue,pink,yellow,pink,orange,purple,red(all colors can be available)
packing	100pcs/bag, 20bags/ctn, 2000pcs/ctn
packing design	All inner box and carton printing can be as your design
carton size	58x25x27cm
Applications	hospital clinical medical personnel and patients. dust-free workshop,, laboratory, food industry, electronic manufactures and so on
MOQ	100,000 PCS

## Pictures of disposable clip cap:





**FAQ:**

**(1) Are you a manufacturer or a trading company?**

We are a manufacturer located in Xiatao city, Hubei. Welcome to visit our factory.

**(2) May I get samples to test quality? Small trail order acceptable?**

Yes, we'd like to send you samples for evaluation. Small trail order is also welcomed.

**(3) Can you offer customized service?**

Yes, customized size and OEM brand available.

**(4) What's your main products?**

Our main products include disposable face mask, disposable bouffant cap, disposable coverall, lab coat, surgical gown, disposable apron, shoe cover, etc.

**(5) What's your delivery time?**

Usually within 15 days for 20feet container.