

## **The property of disposable Customized towels for beauty salon**

**All cotton towel are made of the highest quality materials, providing excellent protection and superior comfort:**

1. Size: 20 x 20 cm each piece or customized
2. Made of cotton material, soft and comfortable, clean and sanitary.
3. Convenient and lightweight , easy to carry out .
4. Can be wet and dry, are great for fragile or sensitive skin
5. Great water absorption, skin-friendly and comfortable.
6. Cotton tissue can instead of face towel and will be much cleaner to dry your face after washing.

## **The Specification of disposable Customized towels for beauty salon**

Name	disposable Customized towels for beauty salon
Material	100% cotton; spunlace nonwoven
Style	Plain, mesh, pearl and customization patterns available
Size	15*20cm, 20*20cm ,or other size can be customized
Color	white
Packing	1roll/portable bag ; 100pcs/pack ,can be customized
Carton Box	All inner bag and carton printing can be as your design
Carton Size	52x38x30cm /52x19x30cm
Applications	personal care, makeup removing, house cleaning, travel use, camping, similar outdoor activities and business occasions.
MOQ	1000 bags
Sample	Can be supplied free for your quality checking within 3 working days
Payment way	T/T, PAYPAL, Westunion, L/C, D/A, D/P, ESCROW
Delivery time	as you quantity
FOB Port	Wuhan/Shanghai/Ningbo/Guangzhou

## **The picture of disposable Customized towels for beauty salon**





COTTON  
ALL CLEAN NATURAL  
CHINA  
100% ORGANIC AND GOOD FEELER  
PROTECTS YOUR SKIN

WET OR DRY

COTTON FACIAL CLEANING TOWEL



# Multi-Purpose



Daily face cleaning



Men's face cleaning



Wipe hands



Make up



No fragrance nor chemicals  
100% biodegradable



Baby's face

