

Description of disposable lab coat:

1. Name: No woven lab coat
2. Application: Hospital, hotel, dustproof place, food/electronic/chemical workshop
3. Material: PP+PE
4. Color: White, green, yellow, blue, pink
5. Style: Elastic/ knitted cuff, round neck/ with collar, front button/ Velcro, Customization is available
6. Size: S-XXXL sizes available
7. Packing: 5PCS/polybag, 50PCS/carton and as you require
8. 1*20" FCL: 50000PCS (tolerance as different sizes and materials)
9. Mini quantity: 10, 000PCS
10. Certificate: CE, FDA, ISO13485

Features of non-woven disposable jackets :

1. Isolation gown is made of excellent non-woven material from biggest supplier
2. Lint-free, waterproof, good tensile strength, comfortable
3. Breathable, excellent to dissipate body heat, against splatter of

Styles of disposable lab coat:

- Collar, open and buttons
- Collar, elastic and buttons
- Collar, open and zipper
- Collar, elastic and zipper
- Collar, open and Velcro
- Collar, elastic and Velcro
- Knitted collar and cuffs, buttons
- Knitted collar and cuffs, zipper

Details picture of disposable dental lab coat:



