

## Description of disposable hospital gown:

1. Material: SMS (Spunbond-Meltblown-Spunbond Non-wovens)
2. SMS barrier fabric is a unique trilaminate construction that offers a high tensile strength and toughness that is also soft, drapable, and easy to work with. SMS provides a fluid and particulate barrier to protect the worker without sacrificing mobility and comfort.
3. Spunbond layer is consist of continuous filament, with good tensile strength and elongation, Melt blown layer is consist of continuous micro fiber
4. Characteristics: Anti-pull, environmentally-friendly, good tensile strength, feeling soft, non-toxic, water-resistant, air- permeable

## The design details for disposable hospital gown :

1. Material: SMS
2. Cuff : elastic or knitted cuff; Neck: ties or Velcro on the back of the neck
3. Waist: ties on the waist; With Paper card or not
4. Color: blue, white, green, yellow, pink etc.
5. Other materials and specifications are also available, can make according to customers' requirements.

## Application:

Can be used in hospital, laboratory, food industry, any where need protection. Prevent and isolate dust, particle, alcohol, blood, bacteria and virus invading. Can be used as a visit coat .Prevent te cross infection in medical treatment at sanitation field. Economical and disposable. Accord with the hygienic standard.

## The Specification of disposable SMS gown:

Material	SMS
Color	White,green,blue,yellow,etc.
Cuff	Knitted cuff or elastic cuff,
Neck	Ties or velcro
Waist	Tie on
Sleeves	Long sleeves/short sleeves/non sleeves
Sterilization	EO Sterilization or non-sterilize
Reinforced layer	as your requirements
Size	S-6XL, or as your request
Minimum order quantity	5,000 pcs
Price:	Negotiable
Payment terms:	T/T, L/C, D/A, D/P, Western Union, ESCROW, Paypal
Supply ability:	1*40 HQ can be finished in 15 days
Delivery time:	According to the quantity, normally, 30 days after confirmed the order

## The size chart of the disposable SMS gown:

Size	S	M	L	XL
Height (CM)	110	115	120	125
chest circumference(CM)	30	137	140	150

Details picture of disposable hospital gown:

