

Product details

Advantage of disposable long sleeve isolation gown:

1. Samples and customization available.
2. Professional QC to control the product quality. All the products comply with international quality standard such as CE,ISO13485, FDA.
3. Have enough experience and gained a high reputation from our customers all The world.Such as America, South America, Europe, Southeast Asia, The Middle East.
4. Quick payment, prompt delivery and competitive price.
5. Stable quality and best after-sale service.
6. OEM service is available.

Specification of disposable long sleeve isolation gown

Name	Doctor White Surgical Disposable Hospital Gowns
Material	Polypropylene,PP Nonwoven Fabric
Style	Non-sterile, standard and reinforced style, standard and ultrasonic, with and without towels. Tie on collar, ties on waist, knit cuffs
Size	S(110*130cm), M(115*137cm), L(120*140cm), XL(125*150cm), 2XL(130*150cm) or any other customized sizes
Grammage	18g
Color	white.
Packaging	non sterile,10 pcs/bag, 50pcs/ctn or 100pcs/ctn,
Container loading	20'FCL:680CTNS,40'HQ:1700CTNS
MOQ	5000 PCS
Delivery	One 40HQ container can be finished within 15-25days
Store	Store in dry, humidity below 80%, ventilated, non-corrosive gases warehouse

Product display



orient Honest
MEDICAL INTERNATIONAL






orient Honest
MEDICAL INTERNATIONAL



Our business process



The regular shipment way



FAQ

Q: 1. Are you a manufacturer or a trading company, and when was your company established?

A: Yes, we are directly manufacturer, our factory was established in 2003, we have our own exported rights company in 2013.

Q: 2. What is the annual turnover for your company?

A: Normally it is around USD18 million-20 million, it is still in the growing.

Q: 3. Do your company have any quality management system such as ISO, EN or CE certificate(s)?

A: Yes, we have CE, ISO13485, FDA, REACH and NELSON report.

Q: 4. Can I get the free samples for quality checking? How long time can I get it?

A: Yes, our samples can be supplied free within 3-7 days.

