

Product details

Features of disposable long sleeve lab coats

- 1.Yellow color,Hand sewing, with kintted cuff and collar,with three pockets.
- 2.Non-irritating,antibacterial property ,good tensile strength,Lint-free,Superior breathability and anti-static,impact resistance,soft,comfortable.

Packing of disposable long sleeve lab coats

Normal packing: 5 pc/bag, 10bags/carton, or 1pc/bag,50pcs/carton.
Your logo can be printed on the inner box and outer carton.

Size Chart

LAB COAT SIZE CHART

(Customized size is welcome)

| CM | S | M | L | XL | XXL | XXXL | XXXXL | XXXXXL |
|--------------|-----|-----|-----|-----|-----|------|-------|--------|
| HEIGHT | 100 | 103 | 107 | 110 | 112 | 113 | 120 | 120 |
| CHEST LENGTH | 132 | 135 | 137 | 140 | 142 | 147 | 150 | 158 |

Specification of disposable long sleeve lab coats

| | |
|--------------------------|---|
| Product: | disposable long sleeve lab coats,disposable surgical coat for hospital,three pocket yellow lab coat |
| Material: | PP non-woven |
| Grammage: | 45gsm |
| Style: | long or short sleeves; elastic or knit cuff ; lapel / stand / round / knitted collar optional with / without pocket; with snap button or Velcro closure |
| Application: | for medical, pharmaceutical, laboratorial, workshop sector |
| Color: | yellow, etc, optional |
| MOQ | 5000pcs |
| Packaging: | 5 or 10pcs/bag; 50 or 100pcs/ctn; |
| Storage: | Stored in dry, humidity below 80%, ventilated, non-corrosive gases warehouse |
| Terms of payment: | T/T, PayPal, Western union, L/C,D/A,D/P,ESCROW, |
| Other service: | OEM and Customized Logo/brand |
| Certification: | CE, ISO13485 |
| Delivery time | Capacity: 10,000pcs/day, 1*40HQ can be finished within 10-15 days; |

Product display



orient Honest
MEDICAL INTERNATIONAL





orient Honest
MEDICAL INTERNATIONAL






orient Honest
MEDICAL INTERNATIONAL



Our business process



The regular shipment way



FAQ

Q: 1.Are you a manufacturer or a trading company, and when was your company established?

A: Yes, we are directly manufacturer, our factory was established in 2003, we have our own exported rights company in 2013.

Q: 2.What is the annual turnover for your company?

A: Normally it is around USD18 million-20 million, it is still in the growing.

Q: 3. Do your company have any quality management system such as ISO, EN or CE certificate(s)?

A: Yes, we have CE,ISO13485,FDA. REACH and NELSON report.

Q: 4.Can I get the free samples for quality checking? How long time can I get it?

A: Yes, our samples can be supplied free within 3-7 days.