

disposable nonwoven surgical caps for doctor

Specification:

| | |
|---------------|--|
| Material | PP, SPP, SBPP, SMS, SMMS, Microporous |
| Color | White, blue, green, yellow, red.....etc |
| Fabric weight | PP 10-30gsm, SMS 30-65gsm, Microporous45-65gsm |
| Type | ties on |
| Elastic | Single/double |
| Size | 18",19",20",21",24"....etc |
| Certificate | CE,ISO13485,FDA |
| Packaging | 100/bag,2000pcs/ctn, or as your requirements |
| Sample | Sample can be finished in 2 days |
| Model | disposable surgical cap |

Advantages:

- 1.Directly manufacturer, disposable mop cap, any colors.
- 2,ISO,CE, good protect, latex free ,easy breath, eco friendly,
- 3.Capacity: 200,000pcs/day, 1*40HQ can be finished within 10-15 days;
4. OEM packing design, service in 24 hours.