

Specification of disposable sterile surgical gown:

Name:	disposable sterile operating gown
Material:	SMS
Grammage:	30-60 gsm
Closure:	with belt closure, or with Velcro closing
Application:	for medical, surgical, pharmaceutical, laboratorial, industrial sector
Size:	S/M/L/XL/XXL/XXXL
Color:	Dark Blue, etc, optional
Packing:	Normal packing:1pcs/polybag,50pcs/carton
Storage:	Stored in dry, humidity below 80%, ventilated, non-corrosive gases warehouse
Certification:	CE, ISO13485

Size Chart of disposable sterile gown :

SURGICAL / ISOLATION GOWN SIZE CHART (Customized size chart is welcome)							
CM	S	M	L	XL	XXL	XXXL	
HEIGHT	115	115	120	125	130	135	
CHEST LENGTH	127	137	140	145	150	155	

Styles of disposable sterile gown:

- A.standard operation gown with knitted cuff
- B.standard operation gown with elastic cuff
- C.half reinforced operation gown with knitted cuff
- D.full reinforced operation gown with knitted cuff

Advantages of disposable sterile gown:

- 1.SMS material: Lint-free, waterproof, good tensile strength, comfortable.
- 2.Light-weight, Non-irritating, antibacterial property ,good tensile strength, Superior breathability and anti-static, impact resistance, soft, comfortable.
- 3.Provide an excellent barrier to potentially infectious blood and other bodily fluids. against splatter of chemical material and solid dust.
- 4.Disposable SMS surgical gown mainly for doctor for operation, also available to use in the isolation room's doctors and visitors.
- 5.Disposable SMS surgical gown is the most conventional medical apparel, can reach the health protection standards in the general case.
- 6.SMS : "Three defenses" (waterproof, anti-alcohol and anti-blood spatter) to comply with domestic and

international medical standards.

Picture of disposable operate gown:

