

Disposable surgical caps Introduction:

Disposable surgical bouffant caps are made of soft PP nonwoven material, sewing and automatic style. Disposable surgical caps nonwoven nurse caps are comfortable for personal protection in hospital or other medical place.

Also the "bouffant" style surgical caps is more comfortable for those with longer hair. This is a good choice because it has more room and doesn't hold sweat against the scalp that can promote irritation and itching. A surgeon can't be worrying about scratching his or her head while performing surgery.

Disposable surgical bouffant caps Advantages:

1. Capability: repels fluids from the outsides.
2. Material: PP nonwoven fabric.
3. Sanitation and quality in accordance with the CE, ISO13485
4. Customers' logos is accepted.
5. Different sizes and colors available.
6. Complete hair covering made from soft comfortable, thermally bonded light polypropylene fibre with an encapsulated elastic edge.
7. Cap is compactly folded in strip easily open to shape.
8. With single or double elastic available, in pleated or clip shape with various colours available.

Disposable surgical caps Indications:

1. Suited to electronic manufactures, restaurants, food processing, school, hospital, etc. Widely used in clinic, hospital restaurant, food shop.
2. Sterile disposable caps are used for hospital clinical medical personnel and patients.
3. Non-sterile disposable caps can be used for dust-free workshop,, hospital, laboratory, food industry, electronic manufactures and so on.

Disposable surgical caps Specifications:

Material	PP non woven
Weight	10-30g
Style	Elastic (latex free), single and double elastic style, with peak
Size	18",19",21",24"28"
Color	blue (all colors can be available)
packing	100pcs/bag, 20bags/ctn, 2000pcs/ctn
packing design	All inner box and carton printing can be as your design
carton size	58x25x27cm
Applications	hospital clinical medical personnel and patients. dust-free workshop,, laboratory, food industry, electronic manufactures and so on
MOQ	100,000 PCS

Details picture of disposable surgical caps:





orient Honest
MEDICAL INTERNATIONAL

