

All [disposable surgical gown](#) are made of the highest quality materials, providing excellent protection and superior comfort:

- Non-linting materials are smooth and comfortable to eliminate skin irritation
- Highly breathable construction layers remain cool, even during long cases
- Various Style can be customized according to your different requirement

Fabric Material	SMS material
Style	With Velcro on collar , tie on waist knitted cuffs
Dimensions size	X-3XL
Fabric weight	30-60g, could be as your request
Color	Blue color
Certificate	ISO13485, FDA
MOQ:	10000 pcs
Packing design	can be as your design
Packing	1pc/plastic paper bag,50pcs/ctn, or as your requirement
Capability	60,000pcs per day
Delivery time	7-30 days
Applications	laboratory, clinic, food industry, dust-free workshop,hospital

The picture of disposable surgical gown





orient Honest
MEDICAL INTERNATIONAL



Ultrasonic seams
More firm and durable



cover main parts with film



Tie-On closure



Hook & Loop
closure



Elastic cuff



Knitted cuff

Multi materials for your choice



PP



SMS



Microporous



