

disposable white lab coats introduction :

1. PP non woven material with anti-static
2. 20-60g/m² weight
3. Elastic/knitted/open cuff
4. long/short sleeve
5. Directly manufacturer, disposable lab coats, any colors
6. Capacity: 10,000pcs/day, 1*40HQ can be finished within 10-15 days
7. Good barrier to fluid and breathable, OEM packing design, service in 24 hours
8. We have CE,FDA and ISO 9001 certification

Disposable white Lab Coat features:

1. Choose of the best non-woven providing reliable and selective barriers to bacteria, blood and other fluids.
2. Lint-free waterproof, good tensile strength, comfortable, breathable. Excellent to dissipate body heat, antistatic, against splatter of liquid
3. Great anti-bacterial effects, protecting workers against exposure to various particulates from pollutants and insecticides etc.

Disposable Lab Coat material features:

PP disposable lab coat made from hydrophobic polypropylene material, Latex-free; abrasion-resistant; low lint; with a high level of fluid repellency.

Disposable visitor Coat Specifications:

Product:	disposable white lab coats
Material:	PP non-woven
Grammage:	20-50gsm
Style:	long or short sleeves; elastic or knit cuff ; lapel / stand / round / knitted collar optional with / without pocket; with snap button or Velcro closure
Application:	for medical, pharmaceutical, laboratorial, workshop sector
Color:	Blue,,White,Dark Blue, etc, optional
MOQ	5000pcs
Packaging:	5 or 10pcs/bag; 50 or 100pcs/ctn;
Storage:	Stored in dry, humidity below 80%, ventilated, non-corrosive gases warehouse
Terms of payment:	T/T, PayPal, Western union, L/C,D/A,D/P,ESCROW,
Other service:	OEM and Customized Logo/brand
Certification:	CE, ISO13485
Delivery time	Capacity: 10,000pcs/day, 1*40HQ can be finished within 10-15 days;

Details picture of disposable white lab coat:

