

## Product details

### Description of SF anti-static uniform lab coat

1. Material:SF50G
2. Color: White,
3. Style: Elastic cuff, with collar, front button,with pockets
4. Size: S-XXXL sizes available
5. Packing: 5PCS/polybag, 50PCS/carton and as you require
- 6.1\*20" FCL: 50000PCS (tolerance as different sizes and materials)
7. Mini order quantity: 5000PCS
8. Certificate: CE, FDA, ISO13485

### Our advantage:

- 1.Directly manufacturer, disposable lab coats, any colors.
- 2,ISO,CE, PP,PP/PE,SMS,SMMS, Microporous,
- 3.Capacity: 10,000pcs/day, 1\*40HQ can be finished within 10-15 days;
- 4.Good barrier to fluid and breathable, OEM packing design, service in 24 hours.
- 5.Disposable lab coat is made from high-performance isolation protective material.

### Specifications of SF anti-static uniform lab coat

Material	Microporous film SF
Grammage	50g ,Could as your request
Style	With collar, Snaps in front,elastic cuff, with pockets
Size	S-6XL
Color	white□ etc.
Packing	5pcs/box,50 bags/ctn, or as your requirement
Packing design	all inner bag and carton printing can be as your design
Applications	dustfree workshop, hospital, laboratory, food industry, electronic manufactures and so on

## Product display



orient Honest  
MEDICAL INTERNATIONAL






orient Honest  
MEDICAL INTERNATIONAL



## Our business process



## The regular shipment way



## FAQ

**Q: 1. Are you a manufacturer or a trading company, and when was your company established?**

A: Yes, we are directly manufacturer, our factory was established in 2003, we have our own exported rights company in 2013.

**Q: 2. What is the annual turnover for your company?**

A: Normally it is around USD18 million-20 million, it is still in the growing.

**Q: 3. Do your company have any quality management system such as ISO, EN or CE certificate(s)?**

A: Yes, we have CE, ISO13485, FDA, REACH and NELSON report.

**Q: 4. Can I get the free samples for quality checking? How long time can I get it?**

A: Yes, our samples can be supplied free within 3-7 days.

