

Product details

Description of SF anti-static uniform lab coat

1. Material:SF50G
2. Color: White,
3. Style: Elastic cuff, with collar, front button,with pockets
4. Size: S-XXXL sizes available
5. Packing: 5PCS/polybag, 50PCS/carton and as you require
- 6.1*20" FCL: 50000PCS (tolerance as different sizes and materials)
7. Mini order quantity: 5000PCS
8. Certificate: CE, FDA, ISO13485

Our advantage:

- 1.Directly manufacturer, disposable lab coats, any colors.
- 2,ISO,CE, PP,PP/PE,SMS,SMMS, Microporous,
- 3.Capacity: 10,000pcs/day, 1*40HQ can be finished within 10-15 days;
- 4.Good barrier to fluid and breathable, OEM packing design, service in 24 hours.
- 5.Disposable lab coat is made from high-performance isolation protective material.

Specifications of SF anti-static uniform lab coat

Material	Microporous film SF
Grammage	50g ,Could as your request
Style	With collar, Snaps in front,elastic cuff, with pockets
Size	S-6XL
Color	white□ etc.
Packing	5pcs/box,50 bags/ctn, or as your requirement
Packing design	all inner bag and carton printing can be as your design
Applications	dustfree workshop, hospital, laboratory, food industry, electronic manufactures and so on

Product display



orient Honest
MEDICAL INTERNATIONAL





orient Honest
MEDICAL INTERNATIONAL



Our business process



The regular shipment way



FAQ

Q: 1. Are you a manufacturer or a trading company, and when was your company established?

A: Yes, we are directly manufacturer, our factory was established in 2003, we have our own exported rights company in 2013.

Q: 2. What is the annual turnover for your company?

A: Normally it is around USD18 million-20 million, it is still in the growing.

Q: 3. Do your company have any quality management system such as ISO, EN or CE certificate(s)?

A: Yes, we have CE, ISO13485, FDA, REACH and NELSON report.

Q: 4. Can I get the free samples for quality checking? How long time can I get it?

A: Yes, our samples can be supplied free within 3-7 days.

