

1. Material:PP non woven/PP+PE/SMS
2. Style :single elastic,double elastic
3. Gram weight:10g-20g
4. Color :White,Blue,Green Etc.
5. Certificate:CE ,ISO 13485,FDA approved

#### Product Parameter

Material	SPP/PP+PE/SMS
Weight	10-20gms
Style	single elastic,double elastic
Size	18',19',20',21'
Color	white, green, blue Etc
packing	100pcs/bag,2000pcs/carton
packing design	All inner bag and carton printing can be as your design
carton size	Standard
Applications	Hospital ,SPA ,Beauty and Etc.
MOQ	100000 PCS
Sample	can be supplied free for your quality checking within one week
Payment term	T/T,L/C
Delivery time	as you quantity(we can finish 100,000 pcs per day)
FOB port	Wuhan/Shanghai/Ningbo/Guangzhou

#### Product Advantage

Supper soft and comfortable

Dust proof

Water resistance

Environmentally friendly Material

#### Product Picture

#### Lead Time

In 15-25 days after confirming receive deposit

Productivity

100,000 pcs per day

Packing

100pcs/bag,2000pcs/carton

100pcs/bag,1000pcs/carton

Shipping Port

Wuhan/Shanghai/Ningbo Port or Any other port in China