

## [headwear cross kf94 mask adults multicolored for adults](#)

**All masks are made of the highest quality materials, providing excellent protection and superior comfort:**

- Non-linting materials are smooth and comfortable to eliminate skin irritation
- Highly breathable construction layers remain cool, even during long cases
- Super soft earloop makes comfortable wearing for long time

## [headwear cross kf94 mask adults multicolored for adults](#)

<b>Fabric Material</b>	Layer 1: PP non woven (customized printing) Layer 2: high filtration paper Layer 3:high filtration paper Layer4: pp non woven
<b>Style</b>	4ply ,with round ear loop
<b>Dimensions size</b>	20*8cm
<b>Fabric weight</b>	30gsm + 25gsm + 25gsm+25gsm, could be as your request
<b>Color</b>	Colorful mask
<b>Certificate</b>	CE, ISO13485, FDA
<b>MOQ:</b>	20000pcs
<b>Packing design</b>	can be as your design
<b>Packing</b>	50 pcs/box,2000 pcs/ctn, or as your requirement
<b>Capability</b>	1 million pcs per day
<b>Delivery time</b>	7-30 days
<b>Applications</b>	hospital, laboratory, clinic, food industry, dust-free workshop

## [The picture of headwear cross kf94 mask adults multicolored for adults](#)



orient Honest  
MEDICAL INTERNATIONAL



orient Honest  
MEDICAL INTERNATIONAL





 orient Honest  
MEDICAL INTERNATIONAL

