

## Description of hospital use nonwoven lab coat :

1. Material: PP50g
2. Color: green,
3. Style: knitted cuff/collar, front button
4. Size: S-XXXL sizes available
5. Packing: 5PCS/polybag, 50PCS/carton and as you require
- 6.1\*20" FCL: 50000PCS (tolerance as different sizes and materials)
7. Mini quantity: 10, 000PCS
8. Certificate: CE, FDA, ISO13485

## Features of hospital use nonwoven lab coat:

1. Hand sewing, green color, knitted cuff and collar with pocket.
2. Non-irritating, antibacterial property, good tensile strength, Lint-free, Superior breathability and anti-static, impact resistance, soft, comfortable.
3. Against splatter of chemical material and solid dust.
4. Usually nonwoven lab coat with the effect of dust-free and anti-virus be widely used public place, medical industry, environment cleaning, electronic factory, food processing, school and hospital.

## LAB COAT SIZE CHART:

Lab coat size chart (Customized size is welcome)								
CM	S	M	L	XL	XXL	XXXL	XXXXL	XXXXXL
HEIGHT	100	103	107	110	112	113	120	120
CHEST LENGTH	132	135	137	140	142	147	150	158

## Specification of hospital use nonwoven lab coat:

**Name:** hospital use nonwoven lab coat, High Quality Medical green Lab Coat  
non woven lab coat manufacturer in China

**Material:** PP non-woven

**Grammage:** 50gsm

<b>Style:</b>	long or short sleeves; elastic or knit cuff ; lapel / stand / round / knitted collar optional with / without pocket; with snap button or Velcro closure
<b>Application:</b>	for medical, pharmaceutical, laboratorial, workshop sector
<b>Color:</b>	green or custom optional
<b>MOQ</b>	5000pcs
<b>Packaging:</b>	5 or 10pcs/bag; 50 or 100pcs/ctn;
<b>Storage:</b>	Stored in dry, humidity below 80%, ventilated, non-corrosive gases warehouse
<b>Terms of payment:</b>	T/T, PayPal, Western union, L/C,D/A,D/P,ESCROW,
<b>Other service:</b>	OEM and Customized Logo/brand
<b>Certification:</b>	CE, ISO13485
<b>Delivery time</b>	Capacity: 10,000pcs/day, 1*40HQ can be finished within 10-15 days;

### **Pictures of hospital use nonwoven lab coat:**



orient Honest  
MEDICAL INTERNATIONAL





orient Honest  
MEDICAL INTERNATIONAL

