

Description of hospital use nonwoven lab coat :

- 1. Material: PP50g
- 2. Color: green,
- 3. Style: knitted cuff/collar, front button
- 4. Size: S-XXXL sizes available
- 5. Packing: 5PCS/polybag, 50PCS/carton and as you require
- 6.1\*20" FCL: 50000PCS (tolerance as different sizes and materials)
- 7. Mini quantity: 10, 000PCS
- 8. Certificate: CE, FDA, ISO13485

Features of hospital use nonwoven lab coat:

- 1.Hand sewing, green color,knitted cuff and collar with pocket.
- 2.Non-irritating,antibacterial property ,good tensile strength,Lint-free,Superior breathability and anti-static,impact resistance,soft,comfortable.
- 3.Against splatter of chemical material and solid dust.
- 4.Usually nonwoven lab coat with the effect of dust-free and anti-virus be widely used public place,medical industry,environment cleaning,electronic factory,food processing,school and hospital.

LAB COAT SIZE CHART:

Lab coat size chart (Customized size is welcome)								
CM	S	M	L	XL	XXL	XXXL	XXXXL	XXXXXL
HEIGHT	100	103	107	110	112	113	120	120
CHEST LENGTH	132	135	137	140	142	147	150	158

Specification of hospital use nonwoven lab coat:

Name:	hospital use nonwoven lab coat, High Quality Medical green Lab Coat□non woven lab coat manufacturer in China
Material:	PP non-woven
Grammage:	50gsm

<b>Style:</b>	long or short sleeves; elastic or knit cuff ; lapel / stand / round / knitted collar optional with / without pocket; with snap button or Velcro closure
<b>Application:</b>	for medical, pharmaceutical, laboratorial, workshop sector
<b>Color:</b>	green or custom optional
<b>MOQ</b>	5000pcs
<b>Packaging:</b>	5 or 10pcs/bag; 50 or 100pcs/ctn;
<b>Storage:</b>	Stored in dry, humidity below 80%, ventilated, non-corrosive gases warehouse
<b>Terms of payment:</b>	T/T, PayPal, Western union, L/C,D/A,D/P,ESCROW,
<b>Other service:</b>	OEM and Customized Logo/brand
<b>Certification:</b>	CE, ISO13485
<b>Delivery time</b>	Capacity: 10,000pcs/day, 1*40HQ can be finished within 10-15 days;

**Pictures of hospital use nonwoven lab coat:**



orient Honest  
MEDICAL INTERNATIONAL





orient Honest  
MEDICAL INTERNATIONAL

