

Disposable gown is made from high-performance isolation protective material. Widely used in medical safety protective, aseptic workshop, protective isolation, mining electronics factory, food factory farm animal husbandry biohazard and so on.

With two kinds processes: sewing and ultrasound thermal bonding.

material features □

#### 1. PP gown

It is made from hydrophobic polypropylene material, Latex-free; abrasion-resistant; low lint; with a high level of fluid repellency.

Color: white, green, blue, yellow, orange and so on

Material weight: 30-65gsm.

#### 2. PP+PE gown

It is made from PP+PE material, Latex-free; abrasion-resistant; completely impervious fluid and alcohol repellent.

Color: white, green, blue, yellow.

Material weight: 40-65gsm.

#### 3. SMS gown (Standard and Reinforced):

It is made from hydrophobic SMS/Spunlace material, Latex-free; abrasion-resistant; low lint; with a high level of fluid repellency; a good barrier to blood, body fluids and pathogens.

With standard and reinforced style. Reinforced style with extra protection reinforcement on sleeve & chest which can be completely impervious fluid and alcohol repellent.

Available color: white, blue.

Material weight: 35-65gsm.

#### 4. Microporous gown

Permeability is strong, anti-static; Better permeability, in a variety of organic solvents, acid and alkali corrosion resistance at the same time, with high resistance to impact mechanical properties of the toughness, texture soft, comfortable not combustion-supporting non-toxic non-irritating, harmless to skin.

Available color: white

Material weight: 45-65gsm.

The Specification of Disposable gown

<b>Fabric Material</b>	PP, PP/PE, SMS, SMMS, CPE
<b>Style</b>	tie on/magic tape on the neck , elastic/knit cuffs,
<b>size</b>	S-5XL
Color	white, blue, green, yellow , pink,etc
Certificate	CE, ISO13485, FDA
MOQ:	5000pcs
Packing design	can be as your design
Packing	10 pcs/box,100 pcs/ctn, or as your requirement
Delivery time	Capacity: 10,000pcs/day, 1*40HQ can be finished within 10-15 days
FOB port	Wuhan/Shanghai/Ningbo/Guangzhou

picture





orient Honest  
MEDICAL INTERNATIONAL





orient Honest  
MEDICAL INTERNATIONAL





