

Specifications of lab coats disposable with pockets:

Name	lab coats disposable with pockets,China plant lab coats for sale , lab coats pink wholesales
Material	PP+PE45G
Style	with three pockets,button on the front, knitted cuff and collar
Size	S M L XL 2XL 3XL 4XL 5XL
Color	pink, or red, grey, black,White, green, blue, yellow, etc.
Packing	Non sterile:5pcs/bag,10 bags/ctn or ; Sterile:1pcs/bag, 50pcs/ctn
Packing design	All inner box and carton printing can be as your design
MOQ	5000 pcs
Sample	Can be supplied free for your quality checking within 3 days
Payment way	T/T, PayPal, Western union, L/C,D/A,D/P,ESCROW,
Delivery time	Capacity: 10,000pcs/day, 1*40HQ can be finished within 10-15 days;
FOB port	Wuhan/Shanghai/Ningbo/Guangzhou/Qingdao/Shenzhen

Lab coats disposable with pockets features:

1. Can be customized with different weight of material and color;
2. Can be single/double/knitted collar, zipper/button/Velcro on the front,
3. Can be elastic/knitted/open cuffs; With one or two or three pockets;
4. Soft, light, non-toxic, durable, eco-friendly, economical, Breathable.
5. Fluid repellent, stain resistant and durable;
6. Prevent and isolate dust, particle, alcohol, blood, bacteria and virus invading.
7. Prevent cross-infection in medical treatment at sanitation field.
8. Can be used as a visit coat.

Sample policy of lab coats disposable with pockets:

1. Sample totally based on your requirements. Sampling time:2-5days.

2. Stock samples sampling time: can send right away.
3. Make mass production after you confirm the samples.
4. Shipping: Samples delivered by DHL, FEDEX, TNT, UPS, EMS; the bulk goods shipping by sea.

Our Advantages:

Professional shipping service if necessary

Full series after sales service system.

Welcome any OEM requirements.

FDA, CE, ISO13485 certification.

Picture of lab coats disposable with pockets:





orient Honest
MEDICAL INTERNATIONAL

