

Introduction of lab coats wholesale for children:

Lab coats wholesale for children made by PP non woven material, widely used in school, hospital and laboratory, providing good protection against dust, particles etc. It's economical for single use.

Specifications of lab coats wholesale for children:

Product:	lab coats wholesale for children
Material:	PP non-woven
Grammage:	40gsm,or grams as your request
Style:	long or short sleeves; elastic or knit cuff ; lapel / stand / round / knitted collar optional with / without pocket; with snap button or Velcro closure
Application:	for medical, pharmaceutical, laboratorial, workshop sector
Color:	Blue,White,Dark Blue, etc, optional
MOQ	5000pcs
Packaging:	1 or 5pcs/bag; 50 or 100pcs/ctn;
Storage:	Stored in dry, humidity below 80%, ventilated, non-corrosive gases warehouse
Terms of payment:	T/T, L/C, D/P,Paypal,etc
Other service:	OEM and Customized Logo/brand
Certification:	CE, ISO13485

Our advantage:

- 1.5, 000 square meters absolutely dustless clean workshop for clothing making
2. We stick to strict quality inspection standard on production ensuring high quality.
3. Different design of products is available on your request.

4. Samples and customization available

Pictures of kids disposable lab coats:





orient Honest
MEDICAL INTERNATIONAL

