

lightweight comfortable [disposable bed cover for hospital /SPA](#)

We are [disposable bed cover maufacturer in China](#) .Made from high-performance isolation protective material. Widely used in hospital for emergency bed and surgical bed . [disposable bed cover](#) is used in hotels and massage bed . To prevent the spread of body fluids and blood splash content and isolate dust ,particle , alcohol, blood , bacteria and virus invading .

Material	PP,PP/PE,SMS,SMMS
Style	Flat, elastic ends, square ends, roll
Size	80x190+20cm,90x200+20cm,160x190+20cm,160x200+20cm,180x200+20cm.....
Color	white,green,blue,yellow,pink etc.
Packing	10pcs/box,10 bags/ctn or 5 pcs/bag,50 pcs/bag
Packing Design	can be as your design
MOQ	5000pcs
Samples	can be supplied free for your quality checking within 3 days
Delivery time	Capacity: 10,000pcs/day, 1*40HQ can be finished within 10-15 days
FOB port	Wuhan/Shanghai/Ningbo/Guangzhou

Picture





orient Honest
MEDICAL INTERNATIONAL







orient Honest
MEDICAL INTERNATIONAL

