

[mask kf94 korean 250g for adults multicolored for adults](#)

All masks are made of the highest quality materials, providing excellent protection and superior comfort:

- Non-linting materials are smooth and comfortable to eliminate skin irritation
- Highly breathable construction layers remain cool, even during long cases
- Super soft earloop makes comfortable wearing for long time

[4ply KF94 face mask black and white color for adults](#)

Fabric Material	Layer 1: PP non woven (customized printing) Layer 2: high filtration paper Layer 3:high filtration paper Layer4: pp non woven
Style	4ply ,with round ear loop
Dimensions size	20*8cm
Fabric weight	30gsm + 25gsm + 25gsm+25gsm, could be as your request
Color	Colorful mask
Certificate	CE, ISO13485, FDA
MOQ:	20000pcs
Packing design	can be as your design
Packing	50 pcs/box,2000 pcs/ctn, or as your requirement
Capability	1 million pcs per day
Delivery time	7-30 days
Applications	hospital, laboratory, clinic, food industry, dust-free workshop

- [The picture of 4ply KF94 face mask](#)



orient Honest
MEDICAL INTERNATIONAL



orient Honest
MEDICAL INTERNATIONAL





