

Product details

Application of disposable surgery gowns for hospital:

1. Surgical gown can be used in hospital, laboratory, food industry, electronic manufactures, anywhere need protection;
2. Surgical gown can be used as a visit protection. Prevent to cross infection in medical treatment at sanitation field. Economical and disposable. Accord with the hygienic standard;

Specification of disposable surgery gowns for hospital:


Name	non-woven surgical gowns workshop, disposable surgery gowns for hospital, China disposable PP surgeon gown
Material	Polypropylene, PP Nonwoven Fabric
Style	Non-sterile, standard and reinforced style, standard and ultrasonic, with and without towels. Tie on collar, ties on waist, knit cuffs
Size	S(110*130cm), M(115*137cm), L(120*140cm), XL(125*150cm), 2XL(130*150cm) or any other customized sizes
Grammage	40g
Color	white.
Packaging	non sterile, 10 pcs/bag, 50pcs/ctn or 100pcs/ctn,
Container loading	20'FCL: 680CTNS, 40'HQ: 1700CTNS
MOQ	5000 PCS
Delivery	One 40HQ container can be finished within 15-25days
Store	Store in dry, humidity below 80%, ventilated, non-corrosive gases warehouse

Product display




orient Honest
MEDICAL INTERNATIONAL





Our business process



The regular shipment way



FAQ

Q: 1. Are you a manufacturer or a trading company, and when was your company established?

A: Yes, we are directly manufacturer, our factory was established in 2003, we have our own exported rights company in 2013.

Q: 2. What is the annual turnover for your company?

A: Normally it is around USD18 million-20 million, it is still in the growing.

Q: 3. Do your company have any quality management system such as ISO, EN or CE certificate(s)?

A: Yes, we have CE, ISO13485, FDA, REACH and NELSON report.

Q: 4. Can I get the free samples for quality checking? How long time can I get it?

A: Yes, our samples can be supplied free within 3-7 days.

The happy time with customers

