

## Specifications:

Item	nonwoven disposable visitor coat□white disposable nonwoven lab coat□PP nonwoven lab coat with buttons
Material	PP non woven fabric
Size	M:115*135CM,L:120*140CM,XL:130*150CM
Weight	17-65 gsm
Color	White, blue, yellow, green, pink and so on
Style	Knitted cuff ,knitted collar ,buttons front
Features	Waterproof, dustproof, anti-static and so on
Application	Hospital,clean room,factory and so on
Package	5 pcs/pe bag,50 bags/ctn or as your requirement
Certificate	CE/ISO13485/FDA

## Description:

- 1.Light and breathable for body.
- 2.Soft hand feeling and comfortable.
- 3.No stimulus to skin.
- 4.Providing reliable barriers to water steam or blood and other fluids.
- 5.It can prevent pollution and the surface of this clothes without stains.
- 6.Materials or other specifications according to customers' requirement.

## Application:

PP non woven disposable lab coat can be used in hospital, laboratory, outdoor antivirus activities, electronic manufactures, any where need protection. Prevent and isolate dust, particle, alcohol, blood, bacteria and virus invading.

Also can be used as a visit protection.Prevent to cross infection in medical treatment at sanitation field. Economical and disposable. According with the hygienic standard.

## Pictures:



orient Honest  
MEDICAL INTERNATIONAL





orient Honest  
MEDICAL INTERNATIONAL

