

Introduction of red color non-woven clip cap:

Red color non-woven clip cap, elastic surgical scrub caps are made of polypropylene material, made by automatic high temperature thermal bonding technology, can be divided into single elastic and double elastic styles.

Product Advantages:

1. Good protective, light weight, latex free.
2. Use the high temperature thermal bonding technology, comfortable to wearing.
3. Any colors can be available.
4. Use eco friendly material, no harm to the human body and natural environment.

Application of elastic surgical scrub caps:

Widely used in food factory, restaurant, and dust-free workshop.

Specifications:

Material	polypropylene, pp nonwoven
Weight	PP 10g
Style	Elastic (latex free), single elastic style
Size	19"
Color	red, pink, blue, white, green, yellow, orange, purple, (all colors can be available)
Packing	100pcs/bag, 20bags/ctn, 2000pcs/ctn
Packing design	All inner box and carton printing can be as your design
Carton size	58x25x27cm
Applications	Hospital clinical medical personnel and patients. Dust-free workshop,, laboratory, food industry, electronic manufactures and so on
MOQ	100,000 PCS

Packing information:

100 pcs/bag, 20 bags/ctn, carton size: 58x25x27cm
1*20 container can load 715ctns, total 1430,000pcs.
1*40 container can load 1737ctns, total 3474,000pcs.

Pictures of non-woven surgical cap:



