

Introduction of disposable PP25g hospital doctor cap with ties on

Doctor cap with ties ons are made of nonwoven fabric, hand made, with two ties on the back of the cap, strong and adjust ties make the cap very comfortable to wear it. Mainly use for hospital doctor operate room etc.

Specifications of disposable PP25g hospital doctor cap with ties on

Material	SBPP/SPP/nonwoven
Weight	PP10-30g,
Style	with ties on the back of the cap
Size	standard size
Color	blue, or as your request
Packing	100pcs/bag, 20bags/ctn, 2000pcs/ctn
Packing design	All inner box and carton printing can be as your design
Carton size	58x25x27cm
Applications	Hospital clinical medical personnel and patients. Dust-free workshop, Laboratory, food industry, electronic manufactures and so on
MOQ	100,000 PCS

Picture of disposable PP25g hospital doctor cap with ties on



orient Honest
MEDICAL INTERNATIONAL





orient Honest
MEDICAL INTERNATIONAL

