

Product details

Details of disposable PE film reinforced surgical gown:

Fabric conforms to **EN13795** standards,

Surgical gown is produced in a **dust free clean room**, which guarantee good quality.

Sterilization could be provide by us or sterile in your local.

ISO13485 and CE certificates are available, contact us for more information!

One packing items:

Description	Quantity	Specification
Surgical gown	1	M-XL
Hand towel	2	30*40cm
Sterilization bag(paper-plastic bag)	1	26*36cm
Wrapping	1	60*60cm

Size chart :

M	120*140cm
L	125*145cm
XL	130*150cm
XXL	135*155cm

Features:

1.SMMS material: Lint-free, waterproof, good tensile strength, comfortable, “Three defenses”

(waterproof, anti-alcohol and anti-blood spatter)

2.Disposable surgical gown mainly for doctor for operation, also available to use in the isolation room’s doctors and visitors.

Picture of disposable PE film reinforced surgical gown:



orient Honest
MEDICAL INTERNATIONAL





Our advantages:

1. 260 workers, 25000 ㎡ dust free workshop,
2. 12 years experience, ISO, CE & FDA certificated,
3. One hour reach Wuhan port,
4. Free samples and quick response,

This is how we keep good quality and reliable cooperation for global valued customers.

Welcome to contact me and call me as the below information:

VIP service
Remey Chi

Mobile: 0086-15623307739

Tel:0086-27 8782 9195

E-mail:sales04@orienthonest.com

Skype:orienthonest04