

Introduction of waterproof isolation gown manufacturer

1. Long sleeves
2. Elastic cuffs.
3. Ties on neck and waist.
4. With sterilized and non-sterilized.
5. PS: different size and color are available.

Specification of waterproof isolation gown manufacturer

Material	PP+PE,etc
Grammage	PP+PE 45g,or grams as your request
Style	Tie on / Velcro on the neck, elastic / cutted cuff, standard or reinforced
Size	120x130cm,130x150cm,150x160cm,150x190cm, 170x170cm
Color	green,etc
Packing	Non sterile,10 pcs/bag,10 bags/ctn. sterile,1 pcs/sterile bag,100 pcs/ctn
Packing design	All inner box and carton printing can be as your design
Carton size	58X28X48 cm
Applications	Hospital clinical medical personnel and patients dust-free workshop, laboratory, food industry, electronic manufactures and so on

Application

Waterproof isolation gown designed to provide a range of protection level depending on the nonwoven material. Made from soft comfortable nonwoven fabric add a pe film material, can safely and environmentally disposed of when soiled requiring no laundering. Help to stop the spread of infections.

Design for hospitals and laboratory environments where the visitor needs to have moderate protection from infectious or chemical contamination.

Sample Policy:

1. Sample according to your design drawing. Sampling time:3-5days .
2. Design your own sample, then remake it.
3. Existing Samples Sampling time: 1day.

4. Make mass production after you confirm the sample.

Quality Control:

1. We will not start to make the products until you confirm the sample.
2. The packing will be water-proof, damp-proof and sealed.

Pictures of waterproof isolation gown manufacturer:





orient Honest
MEDICAL INTERNATIONAL

