

## Product details

### Features of white blue OEM disposable nonwoven lab coat

- 1.Non-irritating,antibacterial property ,good tensile strength,Lint-free,Superior breathability and anti-static,impact resistance,soft,comfortable.
- 2.Against splatter of chemical material and solid dust.
- 3.Usually nonwoven lab coat with the effect of dust-free and anti-virus be widely used public place,medical industry,environment cleaning,electronic factory,food processing,school and hospital.

### Specifications of white blue OEM disposable nonwoven lab coat

Material	PP nonwoven
Fabric weight	40g,or Could as your request
Style	Knitted collar, knitted cuff, without three pockets
Size	S-6XL
Color	white,green,blue,yellow,pink,red, gray,black etc.
Packing	5pcs/box,50 bags/ctn, or as your requirement
Packing design	all inner bag and carton printing can be as your design
Applications	dust free workshop, hospital, laboratory, food industry, electronic manufactures and so on

## Product display



orient Honest  
MEDICAL INTERNATIONAL





orient Honest  
MEDICAL INTERNATIONAL



## Our business process



## The regular shipment way



## FAQ

**Q: 1. Are you a manufacturer or a trading company, and when was your company established?**

A: Yes, we are directly manufacturer, our factory was established in 2003, we have our own exported rights company in 2013.

**Q: 2. What is the annual turnover for your company?**

A: Normally it is around USD18 million-20 million, it is still in the growing.

**Q: 3. Do your company have any quality management system such as ISO, EN or CE certificate(s)?**

A: Yes, we have CE, ISO13485, FDA, REACH and NELSON report.

**Q: 4. Can I get the free samples for quality checking? How long time can I get it?**

A: Yes, our samples can be supplied free within 3-7 days.