

## Product details

### Introduction of yellow pp isolation gown

Yellow pp isolation gown supplier,pp surgical gown for doctor,cheap disposable medical gowns made from hydrophobic polypropylene material, include PP ,PP+PE ,SMS ,microporous ;Latex free; abrasion resistant; low lint; with a high level of fluid repellency.

### Specification of yellow pp isolation gown

Material	PP non woven
Grammage	20g, or 15-65gsm
Color	Yellow ,we can also make white, green, blue, pink ,red, grey, black etc.
Size	120x130cm,130x150cm,150x160cm,150x190cm, 170x170cm
Style	Tie on / Velcro on the neck, elastic / knitted cuff, standard or reinforced
Packing	Non sterile,10 pcs/bag,10 bags/ctn. sterile,1 pcs/sterile bag,100 pcs/ctn
Packing design	All inner box and carton printing can be as your design
Carton size	58X28X48 cm
Applications	Hospital clinical medical personnel and patients. dust-free workshop, laboratory, food industry, electronic manufactures and so on

## Product display



orient Honest  
MEDICAL INTERNATIONAL






orient Honest  
MEDICAL INTERNATIONAL



## Our business process



## The regular shipment way



## FAQ

- 1, Free sample.
- 2, Small order acceptable;
- 3, L/C, T/T, D/P terms available;
4. FOB, CNF, CIF, EXW terms available;
- 5, Competitive and stable price because we are factory.
- 6, Qualified products. 3rd part product inspection is acceptable.
- 7, OEM service. More options on material, pattern, size and packing;
- 8, Timely delivery. The best lead time is 25-30 day for a 40HQ container.
- 9, Professional export service. Documentation, clearance, sea transportation, all in one.
- 10, Good after sale service. We take care complaints as we would like to be your long-term supplier!