

Недвижимость из Disposable White PP lab coat

All cotton towel are made of the highest quality materials, providing excellent protection and superior comfort:

1. Size: S-5XL or customized
2. Made of nonwoven material, breathable and comfortable, clean and sanitary.
3. Convenient and lightweight, easy to carry out.

The Specification of Disposable White PP lab coat

Имя	Disposable White PP lab coat
Material	pp
Стиль	Knitted collar and knit cuffs, with buttons and pockets
Размер	S-5XL
Цвет	синий
packing	5pcs/box, 10 bags/ctn
packing design	all inner box and carton printing can be as your design
Могил	5000 pcs
Образец	can be supplied free for your quality checking within 2 days
Payment way	T/T, PayPal, Western Union, L/C, D/A, D/P, ESCROW,
Срок поставки	Capacity: 10,000pcs/day, 1*40HQ can be finished within 20-25 days;
FOB port	Wuhan/Shanghai/Ningbo/Guangzhou

The picture of Disposable White PP lab coat





