

Indications of disposable shoe covers:

For hospital clinical medical and patients ,dust free work shop ,food factory.

Packing of disposable shoe covers:

100pcs/bag,2000pcs/ ctn. carton size:36*29*54cm.,1*40container can load 240,0000pcs.

| | |
|----------------|------------------------------------------------------------------|
| Name | shoe cover |
| Material | PP,PP+PE,CPE |
| Style | Sewing ,automatic machine made |
| Size | 38x15cm,40x15cm,41x17cm..... |
| Color | White,blue,green,red colors can be available. |
| packing | 100pcs/bag,2000pcs/ctn |
| packing design | all inner box and carton printing can be as your design |
| MOQ | 50000 pcs |
| Sample | can be supplied free for your quality checking within 3 days |
| Payment way | T/T, PayPal, Western union, L/C,D/A,D/P,ESCROW, |
| Delivery time | Capacity: 600,000pcs/day, 1*40HQ can be finished within 30 days; |
| FOB port | Wuhan/Shanghai/Ningbo/Guangzhou |

