

Информация о продукте:

1. Material: Non woven polypropylene 100%
2. Size:80*190+20cm,90*200+20cm,160*190+20cm,160*200+20cm,180*200+20cm or customized size
3. Elastic: 5mm along side width+30cm on length or customized size
4. Colour: Blue; Белый ; other size available
5. Packing: 10pcs/bag, 100pcs/carton

Indications

- 1, Widely used for SPA, Hotel, Massage, Salon, Medical one time use
- 2, To prevent the spread of body fluids and blood splash content, and isolate the dust, particle, alcohol, blood, bacteria and virus invading.

Характеристики:

Материал	Non woven polypropylene 100%
Цвет	White, blue
Стиль	elastics on width
Размер	80*190+20cm,90*200+20cm,160*190+20cm,160*200+20cm,180*200+20cm or customized size
packing	10 pcs/box, 100pcs/ctn
Масса	25-65gsm or customized
Упаковка	PE bag coated, then paper carton package or customized
Характерная черта	Durable, soft, disposable
использование	SPA, Hotel, Massage, Salon, Medicinal, Hospital, Home

Наши сервисы:

1. Provide the highest quality products
Professional QC team strictly control each step of the production until loading.
2. The fastest delivery time
If urgent, 1*40HQ can be finished within 5-7days!
Normally, 1*40HQ can be finished within 10-15 days!
Busy season, 1*40HQ can be finished within 25-30 days!
3. The most reasonable price
Continuously invest in improving our skilled workers, products and processes to minimize the cost of each step!
- 4.The sweet and quick service
Our workers will be always so patiently making a response of your any questions within 5-10 minutes.
- 5.We can also help to source connected products.
Besides all kinds of disposable medical products, we can also help your source the connected products as your demand!





