

Nonwoven Clip Cap, Disposable Cap, Mob Cap, Doctor Cap, Bouffant Cap, Nurse Cap, Polypropylene Cap Disposable PP Clip Cap ,Disposable Bouffant Cap

1. Directly manufacturer, disposable clip caps, any colors.
- 2, ISO, CE, good protect, latex free, easy breath, eco friendly,
3. Capacity: 200, 000PCS/day, 1\*40HQ can be finished within 10-15 days;
4. OEM packing design, service in 24 hours.

Nonwoven Disposable Caps, Disposable Nurse Cap, Surgical Bouffant Cap With Single Elastic

It is made of PP non woven, sewing and automatic style.

#### Advantages

1. Good protective, light weight, latex free.
2. Sewing style, Comfortable wearing.
3. White, blue, green, yellow, pink, orange, red colors can be available.

#### Indications

For hospital clinical medical and patients, dust free work shop, food factory.

#### Clip caps packing

100PCS/bag, 2000PCS/carton. Carton size: 58X25X27cm 1\*20 container can load 715ctns, total 1430, 000PCS. 1\*40 container can load 1737ctns, total 3474, 000PCS

Name	Disposable surgical non woven elastic caps
Material	PP 10-30g
Style	Bouffant cap, clip cap, mob cap. mop cap. doctor cap
Size	18",19",20",21",22",24"
Color	white,green,blue,yellow,pink,red,orange etc.
packing	100 pcs/bag,2000pcs/ctn
packing design	all inner bag and carton printing can be as your design
MOQ	100,000 pcs
Sample	can be supplied free for your quality checking within 3 days
Payment way	T/T, PayPal, Western union, L/C,D/A,D/P,ESCROW,
Delivery time	Capacity: 200,000pcs/day, 1*40HQ can be finished within 10-15 days;
FOB port	Wuhan/Shanghai/Ningbo/Guangzhou

The picture of disposable bouffant cap