

1) 17.5 * 9.5cm, 3ply, 18 + 20 + 20G

2) 17.5 * 9.5cm, 3ply, 25 + 25 + 25G

- 1) 17.5 * 9.5cm, 3ply, 18 + 20 + 20G, 50 pieces / box, 40 boxes / pallet, 17.5 * 9.5cm, 3ply, 18 + 20 + 20G, 50 pieces / box, 40 boxes / pallet
- 2) 17.5 * 9.5cm, 3ply, 25 + 25 + 25G, 50 pieces / box, 40 boxes / pallet
- 3) BFE 99%

3) 17.5 * 9.5cm, 3ply, 25 + 25 + 25G

- 1) 17.5 * 9.5cm, 3ply, 25 + 25 + 25G
- 2) PP + PE + PP
- 3) Burgundy
- 4) 17.5 * 9.5cm, 3ply, 25 + 25 + 25G
- 5) 17.5 * 9.5cm, 3ply, 25 + 25 + 25G

4) 17.5 * 9.5cm, 3ply, 25 + 25 + 25G

- 1) EN14683: 2005
- 2) BFE% > 99%
- 3) IIR, 120mmHg


5) 17.5 * 9.5cm, 3ply, 25 + 25 + 25G

17.5 * 9.5cm, 3ply, 18 + 20 + 20G	50 pieces / box, 40 boxes / pallet	52 * 19 * 42cm	1 * 40 HQ (3,280,000 pieces)
17.5 * 9.5cm, 3ply, 25 + 25 + 25G	50 pieces / box, 40 boxes / pallet	51 * 37 * 30 cm	2,400,000 pieces

6) 17.5 * 9.5cm, 3ply, 25 + 25 + 25G



orient Honest
MEDICAL INTERNATIONAL





orient Honest
MEDICAL INTERNATIONAL

