

PP / PP + PE SMS Microporous

- 1, Model: OH-3000
- 2. Material: PP, PP + PE, SMS, Microporous
- 3. Color: ...
- 4, Weight: 50 gsm / 100 gsm / 150 gsm
- 5, Standards: ISO13485, FDA

...

...

- ...
- ...
- ...

Material	PP, PP + PE, SMS, Microporous
Model	OH-3000
Weight	X-5XL
Weight	30-60 gsm
Color	White, Black, ...
Standards	ISO13485, FDA
MOQ:	5000 gsm
Weight	50 gsm / 100 gsm, 2000 gsm / 150 gsm
Weight	60,000 gsm
Weight	7-30 gsm
Standards	ISO13485, FDA, ...

