

The property of Disposable White PP lab coat

Disposable cotton towel are made of the highest quality materials, providing excellent protection and superior comfort:

1. Size: S-5XL or customized
2. Made of nonwoven material, breathable and comfortable, clean and sanitary.
3. Convenient and lightweight, easy to carry out.

The Specification of Disposable White PP lab coat

Product Name	Disposable White PP lab coat
Material	pp
Features	Knitted collar and knit cuffs, with buttons and pockets
Size	S-5XL
Weight	Lightweight
Packing	5pcs/box, 10 bags/ctn
Packing design	all inner box and carton printing can be as your design
MOQ	5000 pcs
Quality Control	can be supplied free for your quality checking within 2 days
Payment way	T/T, PayPal, Western union, L/C, D/A, D/P, ESCROW,
Delivery time	Capacity: 10,000pcs/day, 1*40HQ can be finished within 20-25 days;
FOB port	Wuhan/Shanghai/Ningbo/Guangzhou

The picture of Disposable White PP lab coat



orient Honest
MEDICAL INTERNATIONAL





orient Honest
MEDICAL INTERNATIONAL





orient Honest
MEDICAL INTERNATIONAL

