

PP Nonwoven Round Bouffant Caps

Product Description:

Product Name	PP Nonwoven Round Bouffant Caps
Options	9G, 10G, 12G, 14G, 18G, 20G, 25G
Material	PP Nonwoven (Polypropylene)
Color	White, Black, Blue
Quantity	100pcs/bag, 1000pcs/ctn
Usage	Medical, Laboratory, Cleanroom, Food Processing
MOQ	100,000 pcs
Standards	1. CE 2. ISO13485, FDA





orient Honest  
MEDICAL INTERNATIONAL





