

## Product details

Product description: SMS nonwoven fabric, SMS nonwoven fabric; abrasion-resistant; repellency

Product name: SMS nonwoven fabric

Product details:

Material: SMS nonwoven fabric, SMS nonwoven fabric

Weight: 50G

Grammage: 50G

Usage: Industrial, Laboratory, Workshop

Application: Laboratory, Workshop

MOQ: 5000pcs

MOQ: 5000pcs

Quantity: 5 pieces 10 pieces / roll; 50 pieces 100pcs / roll;

Material: 80% SMS nonwoven fabric

Payment: T / T, PayPal, L / C, D / A, D / P, etc.

Service: OEM nonwoven fabric / SMS nonwoven fabric

CE, ISO13485

Capacity: 10,000 / roll 1 \* 40HQ container 10-15 days

## Product display



orient Honest  
MEDICAL INTERNATIONAL





orient Honest  
MEDICAL INTERNATIONAL



## Our business process



## The regular shipment way



## FAQ

Q: 1. Are your products made in China?

A: Yes, our products are made in China. We have been producing since 2003 and have a good reputation in the market.

Q: 2. What is the price of your products?

A: The price of our products is USD18 per unit, with a minimum order of 20 units.

Q: 3. What are the standards of your products?

A: Our products meet the standards of CE, ISO13485, and NELSON.

Q: 4. How long does it take to receive the products?

A: It takes 3-7 days to receive the products.