

buoufant

buoufant:	buoufant
Material:	PP PU
Weight:	10-30gsm
Color:	White, Black, Grey, Blue, Green, Red, Yellow, Orange, Purple, Silver, Gold, etc.
Size:	18", 19", 20", 21", 24", etc.
Quantity:	100 pieces / roll 10m / roll 20m
Lead time:	3-5 days for stock, 7-10 days for custom, 15-20 days for new design
Payment:	T/T, L/C, Paypal
MOQ:	100 pieces, 80% for sample
Validity:	2 years
Application:	OEM, Home Textiles, etc.
Certification:	CE, ISO13485

buoufant

- buoufant (10-30GSM) PU
-
-
- buoufant CAP /

80% for sample

buoufant



